

Präoxygenierung mit Hamilton - Vapox

Modus

- NIV-ST

Einstellungen

- AF: 6-8 breaths/min
- P_{insp}: 10 cmH₂O
- PEEP: 5 cmH₂O
- Oxygen: 100%
- ETS (expiratory trigger sensitivity): 50%
- Inspiratory flow trigger: 2 l/min
- T_i: 2 s
- P-ramp: 50 ms

Wenn Patient anpöisch wird Start mit PCV+

Quelle: <https://www.hamilton-medical.com/de/Resource-center/Article-pageknowledge-base/22d10959-2780-46dd-81fa-84b33c84d9dd~.html>

Revision #2

Created 2025-11-16 21:26:27 UTC by Andi

Updated 2025-11-16 21:32:16 UTC by Andi